



ADDITIONAL INFORMATION

Weight	0.05 kg
Dimensions	2.5 × 9.5 × 10 cm
Carton Qty	120
Pricing Unit	Box/1
Size	6cm x 1M
Min Order Qty	6
ARTG ID	159198
GMDN Code	44990
Medical Device Class	Class 1s
UNSPSC Code	42311546
Tax Schedule	Taxable
Item Length	100
Item Width	6
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	60



PRODUCT DESCRIPTION

AeroPlast Sensitive is a comfortable protective dressing range engineered for delicate, reactive, and ageing skin. The hypo-allergenic, latex-free construction supports skin compatibility while providing dependable cover for cuts, grazes, and minor post-procedure sites at home, school, workplaces, and care facilities. A soft, conforming spunlace non-woven backing moves with the body to maintain comfort through daily activities and helps reduce friction on fragile skin.

SupaBond adhesive delivers strong, consistent adhesion designed to resist edge lift on high-motion areas such as fingers, knuckles, and joints. The formulation balances holding power with considerate removal for sensitive skin, supporting change schedules where frequent replacement is preferred for hygiene or comfort. A non-adherent contact layer sits over the wound bed to minimise sticking during removal and to support atraumatic inspection of the site.

The wound pad is soft and highly absorbent to manage exudate efficiently and help limit maceration, supporting a cleaner micro-environment between changes. Sterilisation by ethylene oxide provides confidence for first application direct from the pack. The breathable, air-permeable backing promotes airflow while retaining strength and flexibility, sitting smoothly over contours and hard-to-cover spots for a reliable seal with low profile comfort. The cut-to-size format allows trimming to the required length for custom coverage that maximises fit and reduces waste.

- **Strong yet gentle hold:** SupaBond adhesive maintains secure adhesion while supporting sensitive skin.
- **Skin compatible:** hypo-allergenic and latex-free to help minimise irritation.
- **Low-trauma removal:** non-adherent contact layer supports gentler changes.
- **Absorbent protection:** soft pad manages exudate to help reduce maceration.
- **Breathable and conforming:** air-permeable fabric that flexes for comfortable wear.

This dressing strip uses a **sensitive spunlace non-woven backing** that is hypo-allergenic and latex-free, combining strong adhesion with skin compatibility for delicate or reactive skin.

The 6 cm × 1 m strip can be cut to the exact length or shape required, allowing tailored coverage on fingers, limbs, and contoured areas while helping reduce waste and keeping a compact box for first aid kits.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.