



ADDITIONAL INFORMATION

Weight	0.1 kg
Dimensions	10.4 × 16.4 × 3.9 cm
Carton Qty	96
Pricing Unit	Box/100
Size	Assorted
Min Order Qty	12
ARTG ID	159198
GMDN Code	44990
Medical Device Class	Class 1s
UNSPSC	42311546
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	60

PRODUCT DESCRIPTION

The AEROPLAST Transparent range brings to you a transparent aesthetic film blending with every skin tone and providing your wound with a waterproof barrier in a plaster form with a full range of sizes.

Aeroplast™ Plastic
washproof sheer backing that bends and
conforms to different shapes and forms

Supabond™
Hypo-allergenic adhesive

exclusively engineered,
for low irritation and high adhesion

SuperSorb™ Wound
Exudate Pad
super absorption,
reducing skin maceration

Non-stick Barrier
protects wound with a non-adherent barrier,
minimising trauma on removal

Unique transparent plaster backing

- The AeroPlast Transparent film is designed to bend, stretch and conform to body contours and guarantee freedom of movement. The film is aesthetic and waterproof, and its transparency assists in wound inspection.

SupaSorb™ Wound Pad

- An embossed wound pad provides superior wound management by draining blood or exudates into pockets in the wound pad. The Super absorbent wound pad minimises the number of dressing changes whilst protecting the wound from further injury. SupaSorb™ leaves the wound site clean and comfortable, while the non-adherent film reduces the pain and site trauma when changing the dressing.

SupaBond™ Adhesive

- SupaBond adhesive, exclusively engineered for unrivalled adhesion, provides a secure adhesion across the site, preventing lift up or edge curling.

Comprehensive range of sizes

- There is an AeroPlast Fabric dressing for any application with specialized fingertip and knuckle dressings, large and small patches, standards and extra wide plasters.



AEROPLAST TRANSPARENT ASSORTED DRESSINGS BOX/100

AT3041

Contraindications / Precautions

Change dressing daily or more often according to the condition of the wound. As with all adhesive products apply and remove carefully from sensitive or fragile skin. Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner

Assorted Dressings, 5 Sizes