

AEROPLAST TRANSPARENT X-WIDE STRIP 7.5 X 2.5CM BOX/100

AT401





ADDITIONAL INFORMATION

Weight 0.1 kg

Dimensions $3.3 \times 12 \times 10.2 \text{ cm}$

Carton Qty 36

Pricing Unit Box/100

Size 7.2cm x 1.9cm

Min Order Qty 12

ARTG ID 159198 **GMDN Code** 44990

Medical Device Class 1s

UNSPSC 42311546

Tax Schedule Taxable

Item Length7.2Item Width1.9

Item Measurement Unit cm

Is Medical Device Yes

Is Serialized Item No

Is Dangerous Good No

Shelf Life (Months) 60

PRODUCT DESCRIPTION

AEROPLAST Transparent is a clear, discreet plaster range designed for professional performance and everyday comfort. The transparent film blends with a wide range of skin tones while forming a water-resistant barrier against dirt and moisture. Built to bend, stretch, and conform to body contours, it supports natural movement so wearers can continue tasks at work, school, or home with confidence.

SupaBond™ adhesive delivers a secure, edge-to-edge hold that resists lifting and curling, even on high-motion areas such as fingers and knuckles. The formulation is latex-free and hypoallergenic, offering reliable adhesion while supporting skin comfort and clean removal—key for long shifts and frequent hand washing.

The SupaSorb™ wound pad manages exudate efficiently by drawing fluid into embossed pockets. A non-adherent contact layer helps the dressing release smoothly, minimising disruption to the site during changes and supporting user comfort. The clear backing also



AEROPLAST TRANSPARENT X-WIDE STRIP 7.5 X 2.5CM BOX/100

AT401

assists quick visual checks without removing the dressing.

Available across a comprehensive family of shapes and sizes, AEROPLAST Transparent equips first aid kits, workplace cabinets, and retail shelves with a go-to solution that looks neat and stays put. From hospitality and healthcare to logistics and education, it offers dependable protection with a professional finish.

- Waterproof transparency: blends with skin for discreet protection that resists splashes.
- **SupaBond™ strong hold:** secure adhesion helps prevent edge lift and curling.
- Hypoallergenic & latex-free: skin-friendly comfort and clean removal.
- **SupaSorb™ absorbent pad:** channels fluids to keep the site cleaner.
- Flexible, conforming film: moves with hands and joints for reliable wear.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Standard Strip