

AT30125





ADDITIONAL INFORMATION

Weight 0.02 kg

Dimensions $0.5 \times 9 \times 11 \text{ cm}$

Carton Qty 500

Pricing Unit Env/25

Size 7.2cm x 1.9cm

Min Order Qty 20

ARTG ID 159198 **GMDN Code** 44990

Medical Device Class 1s

UNSPSC 42311546

Tax Schedule Taxable

Item Length 7.2

Item Width 1.9

Item Measurement Unit cm

Is Medical Device Yes

Is Serialized Item No

Is Dangerous Good No

Shelf Life (Months) 60

PRODUCT DESCRIPTION

The AEROPLAST Transparent range brings to you a transparent aesthetic film blending with every skin tone and providing your wound with a waterproof barrier in a plaster form with a full range of sizes.

Unique transparent plaster backing

• The AeroPlast Transparent film is designed to bend, stretch and conform to body contours and guarantee freedom of



AT30125

movement. The film is aesthetic and waterproof, and its transparency assists in wound inspection.

Aeroplast™ Plastic washproof sheer backing that bends and conforms to different shapes and forms

Supabond™ Hypo-allergenic adhesive exclusively engineered, for low irritation and high adhesion

SuperSorb™ Wound Exudate Pad super absorption, reducing skin maceration

Non-stick Barrier protects wound with a non-adherent barrier, minimising trauma on removal

SupaSorb™ Wound Pad

o An embossed wound pad provides superior wound management by draining blood or exudates into pockets in the wound pad. The Super absorbent wound pad minimises the number of dressing changes whilst protecting the wound from further injury. SupaSorb™ leaves the wound site clean and comfortable, while the non-adherent film reduces the pain and site trauma when changing the dressing.



AT30125

SupaBond™ Adhesive
 SupaBond adhesive, exclusively engineered for unrivalled adhesion, provides a secure adhesion across the site, preventing lift up or edge curling.
Comprehensive range of sizes
 There is an AeroPlast Fabric dressing for any application with specialized fingertip and knuckle dressings, large and small patches, standards and extra wide plasters.
Contraindications / Precautions
Change dressing daily or more often according to the condition of the wound. As with all adhesive products apply and remove carefu from sensitive or fragile skin. Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

This product is not a substitute for direct medical attention from a healthcare practitioner



AT30125

Standard Strip