



ADDITIONAL INFORMATION

Weight	0.020833 kg
Dimensions	13.5 × 9 × 5.4 cm
Carton Qty	24
Pricing Unit	Box/12
Size	2.5cm x 5M
Min Order Qty	1
ARTG ID	236937
GMDN Code	34830
Medical Device Class	Class 1
UNSPSC	42311552
Tax Schedule	Taxable
Location	AU-Warehouse
Is Medical Device	Yes
Is Serialized Item	F
Groups of Qty	1
Item Length	500
Item Width	2.5
Item Measurement Unit	cm

PRODUCT DESCRIPTION

AEROPLAST Transparent Tape is a flexible and breathable transparent tape suitable for the fixation of medical devices to the skin in hospital and clinical settings.



AEROTAPE TRANSPARENT PLASTIC TAPE 2.5CM X 5M BOX/12

AT25

Transparent Film
Waterproof and Micro-Perforated
for Ease of Use

Supabond™
Hypo-allergenic adhesive
Exclusively Engineered
for Low Irritation and High Adhesion

FEATURES & BENEFITS

Hypoallergenic

AEROPLAST Transparent Tape has a very low irritant risk, and is suitable for use on fragile and delicate skin.

Ease of Use

AEROPLAST Transparent Tape is micro-perforated for easy tear-off on any direction and easy adjustment to size.

Individually Wrapped

AEROPLAST Transparent Tape comes in boxes of individually wrapped rolls ensuring the tape is kept in optimum condition and ready for use

Transparent

AEROPLAST Transparent Tape is transparent allowing for improved aesthetics and the monitoring the cleanliness of the site without removing the product.

Microperforated Transparent Tape (wrapped rolls)