



ADDITIONAL INFORMATION

Weight	0.020833 kg
Dimensions	13.5 × 9 × 5.4 cm
Carton Qty	24
Pricing Unit	Box/12
Size	2.5cm x 5M
Min Order Qty	1
ARTG ID	236937
GMDN Code	34830
Medical Device Class	Class 1
UNSPSC Code	42311552
Tax Schedule	Taxable
Item Length	500
Item Width	2.5
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	36
Warranty Period (Months)	12



PRODUCT DESCRIPTION

AEROTAPE Transparent Plastic Tape delivers dependable fixation for dressings, tubing and light medical devices in hospitals, clinics and

first aid environments. The clear, breathable plastic film supports skin comfort while allowing quick visual checks of the site without lifting the tape. Designed for routine procedures and post-procedure care, it provides reliable hold with a smooth, low-profile finish that sits neatly on the skin.

The adhesive is gentle for sensitive or fragile areas and is latex free for added peace of mind. AEROTAPE conforms over contours to help maintain a tidy seal through everyday movement, reducing edge lift and snagging on clothing. The micro-perforated construction enables a clean bi-directional tear, so clinicians and carers can customise length and width in moments without scissors—ideal when workflow speed matters or access to tools is limited.

Each roll is individually wrapped to keep it clean until needed, supporting organised storage on trolleys, in kits and on ward shelves. The combination of visibility, comfort and easy handling makes this tape a practical choice for home, work and clinical use. From securing IV hubs and wound dressings to anchoring sensors and light tubing, AEROTAPE brings clarity and control to day-to-day clinical tasks.

Built for everyday efficiency, this clear micro-perforated tape streamlines application and checks: size, place and confirm fixation in seconds, while patients benefit from gentle materials and comfortable wear. Whether used during dressing changes, device stabilisation or routine monitoring, AEROTAPE balances performance and skin friendliness so attention stays on care outcomes.

- **Hypoallergenic and latex free:** gentle on sensitive and fragile skin.
- **Micro-perforated easy tear:** customise length and width without scissors.
- **Transparent, breathable film:** monitor the site while maintaining comfort.
- **Individually wrapped rolls:** clean, ready to use, and convenient for kits.
- **Conforming flexibility:** smooth fit helps secure dressings and devices.

This transparent medical tape range features a clear, micro-perforated plastic film designed to secure dressings and light devices while allowing quick visual checks of the skin and dressing edge.

Supplied as **2.5 cm x 5 m** transparent medical tape in a **box of 12 rolls**, suitable for bulk replenishment and efficient kit restocking while maintaining consistent, gentle fixation.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Microperforated Transparent Tape (wrapped rolls)