

# AEROPLAST WOUND CLOSURE STRIPS 6 X 75MM 3 STRIPS/CARD BOX/50

APWC7563





## ADDITIONAL INFORMATION

Weight	0.077 kg
Dimensions	2.5 × 6.6 × 13.6 cm
Carton Qty	100
Pricing Unit	Box/50
Size	6mm x 75mm
Min Order Qty	1
ARTG ID	159198
GMDN Code	44990
Medical Device Class	Class 1s
UNSPSC	42311546
Tax Schedule	Taxable
Item Length	7.5
ltem Width	0.6
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	36

## **PRODUCT DESCRIPTION**

The AEROPLAST Wound Closure strips offer you a painless alternative for the closure of minor wounds without the need for stitching.



## AEROPLAST WOUND CLOSURE STRIPS 6 X 75MM 3 STRIPS/CARD BOX/50

**APWC7563** 

Coated Non-woven Polyester Backing Resists flex and Movement

Supabond™ Hypo-allergenic adhesive exclusively engineered, for low irritation and high adhesion

#### SupaBond<sup>™</sup> Adhesive

• SupaBond adhesive, exclusively engineered for unrivalled adhesion, provides a secure adhesion across the site, preventing lift up or edge curling.

#### **Comprehensive sizes available**

• AeroPlast wound closures are available in sizes from 38 to 100mm long and from widths of 3mm to 100mm, enough to suit any incision or laceration. They can also be cut to size.

### **Neat and Simple**

• AeroPlast Wound closures offer fast and easy application with no need for an anaesthetic. AeroPlast wound closures reduce scarring and eliminate needle puncture marks.

#### **Contraindications / Precautions**

- Aeroplast should not be used on infected wounds unless under the supervision of a healthcare practitioner.
- As with all adhesive products, apply and remove carefully from sensitive or fragile skin.
- Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

### This product is not a substitute for direct medical attention from a healthcare practitioner.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Wound Closure Strips - cards of 3