



ADDITIONAL INFORMATION

Weight	0.027667 kg
Dimensions	8.2 × 12.5 × 5.8 cm
Carton Qty	30
Pricing Unit	Box/6
Size	5cm x 5M
Min Order Qty	1
ARTG ID	236937
GMDN Code	34830
Medical Device Class	Class 1
UNSPSC	42311552
Tax Schedule	Taxable
Item Length	500
Item Width	5
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	36

PRODUCT DESCRIPTION

AEROTAPE Microporous Paper Tape is a flexible and breathable paper tape suitable for securing dressings, especially to facial and fragile skin. AeroTape Non-woven Microporous Paper Tape is recommended for the securing of dressings and medical equipment on sensitive and fragile skin including facial areas.



AEROTAPE WHITE MICROPOROUS PAPER TAPE

5CM X 5M BOX/6

AP350

AeroFilm
Ultra-thin, Waterproof & Transparent Bacterial Barrier

Supabond™
Hypo-allergenic adhesive
Exclusively Engineered
for Low Irritation and High Adhesion

Hypoallergenic

- AEROTAPE Microporous Paper Tape has a very low irritant risk, and is suitable for use on fragile and delicate skin.

Ease of Use

- AEROTAPE Microporous Paper Tape is available in a ready to use dispenser for easy dispensing and tear off

Individually Wrapped

- AEROTAPE Microporous Paper Tape is available in boxes of individually wrapped rolls ensuring the tape is kept in optimum condition.



AEROTAPE WHITE MICROPOROUS PAPER TAPE

5CM X 5M BOX/6

AP350

Microporous Paper Tape