



ADDITIONAL INFORMATION

Weight	0.3 kg
Dimensions	13.2 × 15.3 × 5.4 cm
Carton Qty	30
Pricing Unit	Box/12
Size	1.25cm x 9.1M
Min Order Qty	1
ARTG ID	236937
GMDN Code	34830
Medical Device Class	Class 1
UNSPSC	42311548
Tax Schedule	Taxable
Item Length	910
Item Width	1.25
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	36

PRODUCT DESCRIPTION

AEROTAPE Microporous Paper Tape is a flexible and breathable paper tape suitable for securing dressings, especially to facial and fragile skin. AeroTape Non-woven Microporous Paper Tape is recommended for the securing of dressings and medical equipment on sensitive and fragile skin including facial areas.



AEROTAPE WHITE MICROPOROUS PAPER TAPE WITH DISPENSER 1.25CM X 9.1M BOX/12

AP312D

AeroFilm
Ultra-thin, Waterproof & Transparent Bacterial Barrier

Supabond™
Hypo-allergenic adhesive
Exclusively Engineered
for Low Irritation and High Adhesion

Hypoallergenic

- AEROTAPE Microporous Paper Tape has a very low irritant risk, and is suitable for use on fragile and delicate skin.

Ease of Use

- AEROTAPE Microporous Paper Tape is available in a ready to use dispenser for easy dispensing and tear off

Microporous Paper Tape with Dispenser



AEROTAPE WHITE MICROPOROUS PAPER TAPE WITH DISPENSER 1.25CM X 9.1M BOX/12

AP312D