



ADDITIONAL INFORMATION

Weight	0.46 kg
Dimensions	14.9 × 20.4 × 5.3 cm
Carton Qty	30
Pricing Unit	Box/12
Size	2.5cm x 9.1M
Min Order Qty	1
ARTG ID	236937
GMDN Code	34830
Medical Device Class	Class 1
UNSPSC	42311548
Tax Schedule	Taxable
Item Length	910
Item Width	2.5
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	9341394000629

PRODUCT DESCRIPTION

AEROTAPE Microporous Paper Tape is a flexible and breathable paper tape suitable for securing dressings, especially to facial and fragile skin. AeroTape Non-woven Microporous Paper Tape is recommended for the securing of dressings and medical equipment on sensitive and fragile skin including facial areas.



AEROTAPE WHITE MICROPOROUS PAPER TAPE WITH DISPENSER 2.5CM X 9.1M BOX/12

AP325D

AeroFilm
Ultra-thin, Waterproof & Transparent Bacterial Barrier

Supabond™
Hypo-allergenic adhesive
Exclusively Engineered
for Low Irritation and High Adhesion

Hypoallergenic

- AEROTAPE Microporous Paper Tape has a very low irritant risk, and is suitable for use on fragile and delicate skin.

Ease of Use

- AEROTAPE Microporous Paper Tape is available in a ready to use dispenser for easy dispensing and tear off

Includes Tape Dispenser.

Microporous Paper Tape with Dispenser



AEROTAPE WHITE MICROPOROUS PAPER TAPE WITH DISPENSER 2.5CM X 9.1M BOX/12

AP325D