



## ADDITIONAL INFORMATION

<b>Weight</b>	0.46 kg
<b>Dimensions</b>	14.9 × 20.4 × 5.3 cm
<b>Carton Qty</b>	30
<b>Pricing Unit</b>	Box/12
<b>Size</b>	2.5cm x 9.1M
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	236937
<b>GMDN Code</b>	34830
<b>Medical Device Class</b>	Class 1
<b>UNSPSC Code</b>	42311548
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	910
<b>Item Width</b>	2.5
<b>Item Measurement Unit</b>	cm
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Shelf Life (Months)</b>	36



## PRODUCT DESCRIPTION

AEROTAPE White Microporous Paper Tape provides dependable, skin-friendly fixation for everyday wound care and clinical use. The non-woven paper construction is microporous to allow airflow, helping skin feel more comfortable during wear while maintaining a reliable hold on dressings and light medical equipment. Designed for delicate applications, it is recommended for use on sensitive or fragile skin, including facial areas and other contoured sites where a gentle touch matters.



# AEROTAPE WHITE MICROPOROUS PAPER TAPE WITH DISPENSER 2.5CM X 9.1M BOX/12

AP325D

This tape is hypoallergenic and latex free to support a wide user base across homes, workplaces, and healthcare environments. The smooth, conformable backing drapes easily over curves so it sits neatly without bulk, assisting quick, tidy applications on hands, limbs, and contour-rich areas. Each roll is individually wrapped to help keep it clean until required, which makes replenishment of first aid kits and treatment rooms straightforward.

Practical handling features help streamline routine care. A simple dispenser core supports controlled unwind, allowing the user to measure and apply the desired length with minimal effort during dressing changes. The balanced adhesive profile is gentle on removal to support wearer comfort while still providing dependable fixation for typical, low-load applications such as securing gauze, pads, or light tubing.

- **Skin-friendly comfort:** hypoallergenic and latex free for sensitive and fragile skin.
- **Breathable support:** microporous, non-woven paper backing allows airflow for comfort.
- **Purpose-built fixation:** recommended for securing dressings and light medical equipment.
- **Hygienic readiness:** individually wrapped rolls stay clean until needed.
- **Easy handling:** dispenser core enables controlled unwind and tidy application.

The **Size: 2.5 cm × 5 m** option offers a versatile width for everyday fixation on hands, limbs and curved areas, providing comfortable coverage while remaining easy to handle on a dispenser core.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Microporous Paper Tape with Dispenser