

ADDITIONAL INFORMATION



Weight	4.94 kg
Dimensions	28 × 17 × 26 cm
Carton Qty	2
Pricing Unit	Bag/100
Size	1180mm x 580mm x 4.5mm
Min Order Qty	1
ARTG ID	236926
GMDN Code	10416
Medical Device Class	Class 1
UNSPSC	42171701
Tax Schedule	Taxable
Item Length	118
Item Width	58
Item Height	0.45
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	60

PRODUCT DESCRIPTION

Hypothermia treatment is a critical part of first aid when dealing with trauma induced shock and where the body has insufficient ability to maintain body temperature.

The Aeroplast Emergency Blanket is a standard foil blanket for inclusion in all first aid kits, to reflect body heat and prevent shock.

Emergency Shock Thermal Blanket



AERORESCUE THERMAL BLANKET 140 X 210CM BAG / 100

ATB130