



ADDITIONAL INFORMATION

Weight	2.2 kg
Dimensions	75 × 35 × 10 cm
Pricing Unit	Each
Min Order Qty	1
ARTG ID	236801
GMDN Code	13818
Medical Device Class	Class 1
Tax Schedule	Taxable
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

- 10 rubberised sturdy grip handles
- The soft stretcher unfolds from a compact storage bag into a large sheet
- Heavy duty, rip resistant PVC vinyl that is both easy to clean and reusable
- Vital in the rapid evacuation of non-ambulance patients from hazardous situations, as it is easy to transport and instantly available
- Foldable and stored in a carry bag for portability.

Product Size: 205 x 70cm
Folding Size: 88 x 24 cm
Weight: 7kg

Load Bearing: 159kg

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.