





## **ADDITIONAL INFORMATION**



Weight 4.52 kg

**Dimensions**  $75 \times 35 \times 10 \text{ cm}$ 

Pricing Unit Each

Min Order Qty 1

**Tax Schedule** Taxable

Is Serialized Item No
Is Dangerous Good No



## **PRODUCT DESCRIPTION**

- 10 rubberised sturdy grip handles
- The soft stretcher unfolds from a compact storage bag into a large sheet
- Heavy duty, rip resistant PVC vinyl that is both easy to clean and reusable
- Vital in the rapid evacuation of non-ambulance patients from hazardous situations, as it is easy to transport and instantly available
- Foldable and stored in a carry bag for portability.

Product Size: 205 x 70cm

Folding Size: 88 x 24 cm



## **AERORESCUE SOFT STRETCHER**

STR-S

Weight: 7kg

Load Bearing: 159kg