







ADDITIONAL INFORMATION

Weight 4.52 kg

Dimensions $75 \times 35 \times 10 \text{ cm}$

No

Pricing Unit Each

Min Order Qty 1

Tax Schedule Taxable

Is Serialized Item No

Is Medical Device

Is Dangerous Good No



PRODUCT DESCRIPTION

- 10 rubberised sturdy grip handles
- The soft stretcher unfolds from a compact storage bag into a large sheet
- Heavy duty, rip resistant PVC vinyl that is both easy to clean and reusable
- Vital in the rapid evacuation of non-ambulance patients from hazardous situations, as it is easy to transport and instantly available
- Foldable and stored in a carry bag for portability.

Product Size: 205 x 70cm Folding Size: 88 x 24 cm

Weight: 7kg



AERORESCUE SOFT STRETCHER

STR-S

Load Bearing: 159kg

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.