

AEROSCAR SILICONE SCAR PATCH 7.4 X 4.4CM BOX/10

ASGS10







ADDITIONAL INFORMATION

Weight 0.06 kg

Dimensions $14 \times 9.6 \times 2.4 \text{ cm}$

Carton Qty 20

Pricing Unit Box/10

Size 7.4 x 4.4cm

Min Order Qty 1

ARTG ID 406607
GMDN Code 44898
Medical Device Class 1

Tax Schedule Taxable

Item Length 7.4

Item Width 4.4

Item Measurement Unit cm

Is Medical Device Yes

3 Wedical Device 103

Is Serialized Item No

Is Dangerous Good No

PRODUCT DESCRIPTION

Our AEROSCAR Silicone Scar Patch is a medical-grade silicone patch designed for the temporary management of both existing and new hypertrophic or keloid scars. It can also be used as a prophylactic treatment on closed wounds to prevent hypertrophic or keloid scarring. The patch is made from a soft, flexible silicone material that is gentle on the skin and can be cut to size to fit individual scars. **Important:** Do not use on open wounds. **Benefits of using AEROSCAR Silicone Scar Patch include:**

- Reduces the appearance of hypertrophic and keloid scars
- Prevents scarring on closed wounds
- Soft, flexible silicone material for comfortable wear
- Easy to apply and remove
- Can be cut to size to fit individual scars

Packaged in boxes of 10 units, each measuring 7.4 cm x 4.4 cm, our AEROSCAR Silicone Scar Patches are easy to store and transport, making them perfect for clinics and medical professionals.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the



AEROSCAR SILICONE SCAR PATCH 7.4 X 4.4CM BOX/10

ASGS10

directions for use.