



## ADDITIONAL INFORMATION

<b>Weight</b>	0.0102 kg
<b>Dimensions</b>	9.2 × 7.2 × 2.5 cm
<b>Carton Qty</b>	600
<b>Pricing Unit</b>	Unit/1
<b>Size</b>	One Size
<b>Min Order Qty</b>	100
<b>ARTG ID</b>	236803
<b>GMDN Code</b>	36066
<b>Medical Device Class</b>	Class 1
<b>UNSPSC</b>	42272303
<b>Tax Schedule</b>	Non-Taxable
<b>Item Length</b>	2
<b>Item Width</b>	6.5
<b>Item Height</b>	7.5
<b>Item Measurement Unit</b>	cm
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Shelf Life (Months)</b>	60



## PRODUCT DESCRIPTION

AEROSHIELD and AEROMASK provide a selection of protective and airway management products for pre-hospital treatment on injured patients, and to assist in the safe application of essential CPR.

### CPR Face Shield

A disposable faces shield made with a tough vinyl and a one way valve for safe application of essential CPR.

Disposable Face Shield



# AEROSHIELD DISPOSABLE FACE SHIELD WITH ONE-WAY VALVE

AFS001