



## ADDITIONAL INFORMATION

<b>Weight</b>	0.0102 kg
<b>Dimensions</b>	9.2 × 7.2 × 2.5 cm
<b>Carton Qty</b>	600
<b>Pricing Unit</b>	Unit/1
<b>Size</b>	One Size
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	236803
<b>GMDN Code</b>	36066
<b>Medical Device Class</b>	Class 1
<b>UNSPSC Code</b>	42272303
<b>Tax Schedule</b>	Non-Taxable
<b>Item Length</b>	2
<b>Item Width</b>	6.5
<b>Item Height</b>	7.5
<b>Item Measurement Unit</b>	cm
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Shelf Life (Months)</b>	60
<b>Warranty Period (Months)</b>	12



AEROSHIELD Disposable Face Shield is designed for safe, effective CPR in public, workplace, and training environments. The clear, tough vinyl barrier covers the nose and mouth while allowing visibility for correct positioning. Its practical build supports reliable ventilations without direct mouth-to-mouth contact—an essential addition to professional kits and community first aid resources.

A built-in one-way valve helps direct airflow to the patient while reducing backflow, supporting hygienic rescue breaths when seconds matter. The single-use format aligns with infection control best practice and simplifies post-event clean-up, making replenishment straightforward for organisations that value readiness and compliance.

The latex-free, hypoallergenic construction supports broad user compatibility across workplaces and public settings. Lightweight and compact, it fits easily in first aid kits, vehicles, or responder pouches so help is at hand during emergencies or drills. It is also suitable for manikin practice, helping learners build confidence with equipment that mirrors real use.

Whether equipping a workplace, stocking community responders, or supporting training programs, this dependable CPR barrier delivers the right balance of protection, clarity, and simplicity for rapid response.

- **One-way valve protection:** promotes hygienic rescue breaths and reduces backflow.
- **Durable vinyl film:** clear, tough barrier that conforms for accurate placement.
- **Latex-free and hypoallergenic:** suitable for diverse users and settings.
- **Single-use format:** supports infection control and easy replenishment.

**Type: Face Shield** delivers a transparent barrier with an integrated one-way valve to support hygienic rescue breaths and confident CPR performance in workplaces, community settings, and training.

**Face Shield** with integrated one-way valve for hygienic rescue breaths. This single-use, clear vinyl barrier is supplied as an individual unit, ideal for first aid kits, vehicles, and responder pouches.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Disposable Face Shield