

AEROSHIELD DISPOSABLE FACE SHIELD IN SACHET BOX/10

AFS003



ADDITIONAL INFORMATION

Weight 0.15 kg

Dimensions $6 \times 7.5 \times 10.4 \text{ cm}$

Carton Qty 144

Pricing Unit Box/10

Size One Size

Min Order Qty 1

 ARTG ID
 236803

 GMDN Code
 36066

Medical Device Class 1

UNSPSC 42272303

Tax Schedule Non-Taxable

Is Medical Device Yes
Is Serialized Item No
Is Dangerous Good No

PRODUCT DESCRIPTION

AEROSHIELD and AEROMASK provide a selection of protective and airway management products for pre-hospital treatment on injured patients, and to assist in the safe application of essential CPR.

CPR Face Shield

A disposable face shield made with a tough vinyl and heat welded antimicrobial filter for safe application of essential CPR.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Premium Disposable Face Shield in Sachet



AEROSHIELD DISPOSABLE FACE SHIELD IN SACHET BOX/10

AFS003