



ADDITIONAL INFORMATION

Weight	0.15 kg
Dimensions	6 × 7.5 × 10.4 cm
Carton Qty	144
Pricing Unit	Box/10
Size	One Size
Min Order Qty	1
ARTG ID	236803
GMDN Code	36066
Medical Device Class	Class 1
UNSPSC	42272303
Tax Schedule	Non-Taxable
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	9341394007352



PRODUCT DESCRIPTION

AEROSHIELD and AEROMASK provide a selection of protective and airway management products for pre-hospital treatment on injured patients, and to assist in the safe application of essential CPR.

CPR Face Shield

A disposable face shield made with a tough vinyl and heat welded antimicrobial filter for safe application of essential CPR.

Premium Disposable Face Shield in Sachet



AEROSHIELD DISPOSABLE FACE SHIELD IN SACHET BOX/10

AFS003