

## AEROSHIELD LARGE NITRILE FINGER COTS BOX/50

ASN050-L





## **ADDITIONAL INFORMATION**

Weight 0.38 kg

**Dimensions**  $8.9 \times 8 \times 2.9 \text{ cm}$ 

Carton Qty 216

Pricing Unit Box/50

Size Large

Min Order Qty 1

**ARTG ID** 238354 **GMDN Code** 11719

Medical Device Class 1

**UNSPSC** 42132202

Tax Schedule Taxable

**Is Medical Device** Yes

**Is Serialized Item** No

**Is Dangerous Good** No

## **PRODUCT DESCRIPTION**

AeroShield Nitrile Finger Cots provide dependable fingertip protection for clinical, workplace and first aid environments where hygiene and precision matter. Made from durable nitrile, they combine flexibility with puncture resistance to create a comfortable barrier that supports extended wear. An ambidextrous, anatomically shaped design helps each cot fit naturally on either hand, streamlining application during busy shifts and helping maintain dexterity for delicate tasks.

The powder-free surface delivers a confident wet and dry grip, aiding accurate handling of instruments, components and dressings. The high-visibility blue colour makes damage easier to spot, supporting quick replacement when required. Anti-static properties broaden use to electronics handling and other static-sensitive settings. Each cot is pre-rolled to speed application when time is critical, from first aid response to routine medical examinations, laboratory work and food preparation.

Engineered for reliability and comfort, AeroShield Nitrile Finger Cots are a practical addition to professional kits and workplace safety programs, helping maintain clean touchpoints while allowing precise control. Their consistent performance supports confident use across healthcare, light industrial, hospitality and education settings.

- **Powder-free wet/dry grip:** secure handling and procedural control.
- Durable nitrile barrier: flexible protection with puncture resistance.
- Blue for visibility: easier split detection during use.
- Pre-rolled and ambidextrous: quick, comfortable application on either hand.



## AEROSHIELD LARGE NITRILE FINGER COTS BOX/50

ASN050-L

• Anti-static: suitable for electronics and static-sensitive environments.

**Size: Large** provides a comfortable, secure fit for broader fingers while maintaining dexterity. The nitrile construction offers flexible, puncture-resistant protection with powder-free wet/dry grip for confident handling.

Supplied in a **box of 50**, ideal for regular replenishment in clinics, labs and workstations. Large blue nitrile cots are powder-free, anti-static and pre-rolled, delivering fast application, reliable grip and easy split detection.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Large Finger Stalls Box/50