

AEROSHIELD SMALL NITRILE FINGER COTS BOX/100

ASN100-S





ADDITIONAL INFORMATION

Weight 0.0585 kg

Dimensions $4.2 \times 8 \times 8.1 \text{ cm}$

Carton Qty 100

Pricing Unit Box/100

Size Small

Min Order Qty 1

ARTG ID 238354 **GMDN Code** 11719

Medical Device Class 1

UNSPSC 42132202

Tax Schedule Taxable

Is Medical Device Yes
Is Serialized Item No
Is Dangerous Good No

PRODUCT DESCRIPTION

AeroShield Nitrile Finger Cots provide dependable fingertip protection for clinical, workplace and first aid environments where hygiene and precision matter. Made from durable nitrile, they combine flexibility with puncture resistance to create a comfortable barrier that supports extended wear. An ambidextrous, anatomically shaped design helps each cot fit naturally on either hand, streamlining application during busy shifts and helping maintain dexterity for delicate tasks.

The powder-free surface delivers a confident wet and dry grip, aiding accurate handling of instruments, components and dressings. The high-visibility blue colour makes damage easier to spot, supporting quick replacement when required. Anti-static properties broaden use to electronics handling and other static-sensitive settings. Each cot is pre-rolled to speed application when time is critical, from first aid response to routine medical examinations, laboratory work and food preparation.

Engineered for reliability and comfort, AeroShield Nitrile Finger Cots are a practical addition to professional kits and workplace safety programs, helping maintain clean touchpoints while allowing precise control. Their consistent performance supports confident use across healthcare, light industrial, hospitality and education settings.

- **Powder-free wet/dry grip:** secure handling and procedural control.
- Durable nitrile barrier: flexible protection with puncture resistance.
- Blue for visibility: easier split detection during use.
- Pre-rolled and ambidextrous: quick, comfortable application on either hand.



AEROSHIELD SMALL NITRILE FINGER COTS BOX/100

ASN100-S

• Anti-static: suitable for electronics and static-sensitive environments.

Size: Small provides a snug, comfortable fit for slimmer fingers, preserving dexterity for delicate tasks. The flexible nitrile film and powder-free grip support precise handling while the blue colour aids quick split detection.

Supplied in a **box of 100** to support frequent changeovers and bulk replenishment. The **small size** offers a snug, comfortable fit for slimmer fingers, helping maintain precise control across clinical, lab and food prep tasks.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Nitrile Finger Cots