

AEROSPLINT FINGER ALUMINIUM FOAM SPLINT 9.5 X 4.5CM

AR1014



ADDITIONAL INFORMATION

Weight 0.005 kg

Dimensions $11 \times 6 \times 0.5 \text{ cm}$

Carton Qty 500

Pricing Unit Each

Size 9.5 x 4.5cm

Min Order Qty 1

ARTG ID 236793
GMDN Code 13681
Medical Device Class 1

UNSPSC 42241507

Tax Schedule Taxable

Item Length 9.5

Item Width 4.5

Item Measurement Unit cm

Is Medical Device Yes

Is Serialized Item No

Is Dangerous Good No

PRODUCT DESCRIPTION

Light Weight Flat Splint

Immobilisation is a critical item of care when responding to a casualty to limit the occurrence of further injury. Quality immobilisation using suitable products is integral to a full recovery of skeletal injuries.

Aluminium Foam Splints

An aluminium foam splint is a highly adaptive splint that can be rolled folded or formed into almost any application with the use of a quality support bandage. Available in a rolled or flat form, shrink-wrapped for protection.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the



AEROSPLINT FINGER ALUMINIUM FOAM SPLINT 9.5 X 4.5CM

AR1014

- 1	٠					_	
~	1	rΔ	cti	\sim	nc	tor	use.
u	ı	ıc	υu	v	ııs	101	use.

Aluminium Foam Splint Rolled