



ADDITIONAL INFORMATION

| | |
|-------------------------------|--------------------|
| Weight | 0.2 kg |
| Dimensions | 14 × 2.3 × 23.3 cm |
| Carton Qty | 50 |
| Pricing Unit | Each |
| Size | 91 x 14cm |
| Min Order Qty | 1 |
| ARTG ID | 236793 |
| GMDN Code | 13681 |
| Medical Device Class | Class 1 |
| UNSPSC | 42241507 |
| Tax Schedule | Taxable |
| Item Length | 91 |
| Item Width | 14 |
| Item Measurement Unit | cm |
| Is Serialized Item | No |
| Is Dangerous Good | No |
| Barcode (GTIN/EAN/UPC) | 9341394026322 |

PRODUCT DESCRIPTION

AeroSplint Folded Extra Large Aluminium Foam Splint

Quality immobilisation is critical for full recovery of skeletal injuries.

The AeroSplint is a highly adaptive splint that can be folded to support a fractured or injured limb.

AeroSplint can be folded to better suit the application site, ensuring a secure fit for splinting.

Shrink wrapped for protection the Extra Large AeroSplint offers greater support for larger limbs.

Other sizes are available.



AEROSPLINT FOLDED EXTRA LARGE ALUMINIUM FOAM SPLINT 91 X 14CM

AR1091XL

Aluminium Foam Splint Rolled