



ADDITIONAL INFORMATION

| | |
|-----------------------------|-------------------|
| Weight | 0.023 kg |
| Dimensions | 14.7 × 7.6 × 2 cm |
| Carton Qty | 585 |
| Pricing Unit | Box |
| Min Order Qty | 1 |
| ARTG ID | 426930 |
| GMDN Code | 36295 |
| Medical Device Class | Class 1 |
| Tax Schedule | Taxable |
| Is Medical Device | Yes |
| Is Serialized Item | No |
| Is Dangerous Good | No |
| Shelf Life (Months) | 60 |



PRODUCT DESCRIPTION

AeroStep™ Gel Sore Spots are made from soft gel with a self-adhesive backing. It provides pressure and friction relief and is difficult to see when in the shoe.

DIRECTIONS: Remove film backing. Apply adhesive side where required. Wash with warm water and let dry. Use for more than 16 hours per day is not recommended.

TOPICAL USE ONLY - DO NOT PUT IN CONTACT WITH OPEN WOUNDS. DO NOT CLEAN WITH ALCOHOL. KEEP OUT OF REACH OF CHILDREN.

THIS PRODUCT IS NOT A SUBSTITUTE FOR PROFESSIONAL ATTENTION FROM A HEALTHCARE PRACTITIONER. DIABETICS AND PEOPLE WITH BLOOD CIRCULATION PROBLEMS SHOULD CONSULT WITH A PHYSICIAN BEFORE USE.



AEROSTEP GEL SORE SPOTS 6PK

ASTSS06

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.