



**ADDITIONAL INFORMATION**

<b>Weight</b>	0.02 kg
<b>Dimensions</b>	14.7 × 7.6 × 2 cm
<b>Carton Qty</b>	585
<b>Pricing Unit</b>	Box
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	426930
<b>GMDN Code</b>	36295
<b>Medical Device Class</b>	Class 1
<b>Tax Schedule</b>	Taxable
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Shelf Life (Months)</b>	60



**PRODUCT DESCRIPTION**

AeroStep™ Gel Toe Cap is an elastic, soft gel cover that protects the toes and nails against toe and footwear chafing. As it is moisturising, elastic and soft, it also helps relieve pain from corns and abrasions.

**DIRECTIONS:** Clean and dry area before use. Apply the toe cap on the affected area. Each wear should not exceed 16 hours. Before reuse, wash cap with warm water and let dry.

TOPICAL USE ONLY - DO NOT PUT IN CONTACT WITH OPEN WOUNDS. DO NOT CLEAN WITH ALCOHOL. KEEP OUT OF REACH OF CHILDREN.

THIS PRODUCT IS NOT A SUBSTITUTE FOR PROFESSIONAL ATTENTION FROM A HEALTHCARE PRACTITIONER.

DIABETICS AND PEOPLE WITH BLOOD CIRCULATION PROBLEMS SHOULD CONSULT WITH A PHYSICIAN BEFORE USE.



# AEROSTEP GEL TOE CAP 2PK

ASTTC02

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.