



ADDITIONAL INFORMATION

Weight	0.024 kg
Dimensions	14.7 × 7.6 × 2 cm
Carton Qty	585
Pricing Unit	Box/2
Min Order Qty	1
ARTG ID	427138
GMDN Code	13544
Medical Device Class	Class 1
Tax Schedule	Taxable
Location	AU-Cameron Dr
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	60



PRODUCT DESCRIPTION

AeroStep™ Gel Toe Spreader and Bunion Protector relieves and dampens pressure and friction in the bunion area, while helping correct positioning of the big toe. Made of an elastic and soft, moisturising gel, it also aids the disappearance of corns.

DIRECTIONS: Clean and dry area before use. Apply the pure gel toe tube on the affected area. Each wear should not exceed 16 hours. Before reuse, wash tube with warm water and let dry.

TOPICAL USE ONLY - DO NOT PUT IN CONTACT WITH OPEN WOUNDS. DO NOT CLEAN WITH ALCOHOL. KEEP OUT OF REACH OF CHILDREN.

THIS PRODUCT IS NOT A SUBSTITUTE FOR PROFESSIONAL ATTENTION FROM A HEALTHCARE PRACTITIONER.

DIABETICS AND PEOPLE WITH BLOOD CIRCULATION PROBLEMS SHOULD CONSULT WITH A PHYSICIAN BEFORE USE.



AEROSTEP GEL TOE SPREADER 22 X 35MM BOX/2

ASTTS02

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.