

## AEROSTEP GEL TOE SPREADER AND BUNION PROTECTOR 1PK

ASTBP01



## ADDITIONAL INFORMATION

Weight	0.025 kg
Dimensions	14.7 × 7.6 × 0.2 cm
Carton Qty	585
Pricing Unit	Box
Min Order Qty	1
ARTG ID	426930
GMDN Code	36295
Medical Device Class	Class 1
Tax Schedule	Taxable
Location	AU-Cameron Dr
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	60

## **PRODUCT DESCRIPTION**

AeroStep<sup>™</sup> Gel Toe Spreader separates the first and second toes, holding overlapping toes in their correct position. Made of an elastic, soft, moisturising gel,

it aids in the prevention of bunions.

**DIRECTIONS:** Place between the toes. Remove overnight. Wash with warm water and let dry. Use for more than 16 hours per day is not recommended.

TOPICAL USE ONLY - DO NOT PUT IN CONTACT WITH OPEN WOUNDS. DO NOT CLEAN WITH ALCOHOL. KEEP OUT OF REACH OF CHILDREN.

THIS PRODUCT IS NOT A SUBSTITUTE FOR PROFESSIONAL ATTENTION FROM A HEALTHCARE PRACTITIONER.



## AEROSTEP GEL TOE SPREADER AND BUNION PROTECTOR 1PK

ASTBP01

DIABETICS AND PEOPLE WITH BLOOD CIRCULATION PROBLEMS SHOULD CONSULT WITH A PHYSICIAN BEFORE USE.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.