



## ADDITIONAL INFORMATION

<b>Weight</b>	0.023 kg
<b>Dimensions</b>	14.7 × 70.6 × 2 cm
<b>Carton Qty</b>	585
<b>Pricing Unit</b>	Box/4
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	428730
<b>GMDN Code</b>	61945
<b>Medical Device Class</b>	Class 1
<b>Tax Schedule</b>	Taxable
<b>Location</b>	AU-Cameron Dr
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Shelf Life (Months)</b>	60



## PRODUCT DESCRIPTION

AeroStep™ Multiuse Fleecy Cotton Pad is a self-adhesive sheet that can be cut to different sizes depending on need. Ideal for protection on painful areas due to friction and corns. Used for pressure relief and to make metatarsal protectors for corns, hard skin on the sole, as well as interdigital.

**DIRECTIONS:** Clean and dry application area before use. Cut the pad to the size and shape needed. Remove residue with warm water if needed.

SINGLE USE ONLY. TOPICAL USE ONLY - DO NOT PUT IN CONTACT WITH OPEN WOUNDS. DO NOT CLEAN WITH ALCOHOL. KEEP OUT OF REACH OF CHILDREN.

THIS PRODUCT IS NOT A SUBSTITUTE FOR PROFESSIONAL ATTENTION FROM A HEALTHCARE PRACTITIONER.

DIABETICS AND PEOPLE WITH BLOOD CIRCULATION PROBLEMS SHOULD CONSULT WITH A PHYSICIAN BEFORE USE.



# AEROSTEP MULTIUSE FLEECY COTTON PAD 68 X 100MM BOX/4

ASTCP04

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.