

AEROSTEP MULTIUSE FLEECY COTTON PAD 68 X 100MM BOX/4

ASTCP04





ADDITIONAL INFORMATION

Weight 0.023 kg

Dimensions $14.7 \times 70.6 \times 2 \text{ cm}$

Class 1

Carton Qty 585

Pricing Unit Box/4

Min Order Qty 1

Medical Device Class

ARTG ID 428730

GMDN Code 61945

Tax Schedule Taxable

Is Serialized ItemNoIs Dangerous GoodNo

Shelf Life (Months) 60

PRODUCT DESCRIPTION

AeroStep™ Multiuse Fleecy Cotton Pad is a self-adhesive sheet that can be cut to different sizes depending on need. Ideal for protection on painful areas due to friction and corns. Used for pressure relief and to make metatarsal protectors for corns, hard skin on the sole, as well as interdigital.

DIRECTIONS: Clean and dry application area before use. Cut the pad to the size and shape needed. Remove residue with warm water if needed.

SINGLE USE ONLY. TOPICAL USE ONLY - DO NOT PUT IN CONTACT WITH OPEN WOUNDS. DO NOT CLEAN WITH ALCOHOL. KEEP OUT OF REACH OF CHILDREN.

THIS PRODUCT IS NOT A SUBSTITUTE FOR PROFESSIONAL ATTENTION FROM A HEALTHCARE PRACTITIONER.

DIABETICS AND PEOPLE WITH BLOOD CIRCULATION PROBLEMS SHOULD CONSULT WITH A PHYSICIAN BEFORE USE.



AEROSTEP MULTIUSE FLEECY COTTON PAD 68 X 100MM BOX/4

ASTCP04

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.