

## AEROSTEP PROTECTIVE BUNION SLEEVE 8.5 X 9CM BOX/1

ASTBS01



## **ADDITIONAL INFORMATION**

Weight 0.022 kg

**Dimensions**  $14.7 \times 7.6 \times 2 \text{ cm}$ 

Carton Qty 585

Pricing Unit Box/1

Min Order Qty 1

ARTG ID 426930
GMDN Code 36295
Medical Device Class 1

**Tax Schedule** Taxable

AU-Cameron Dr

Is Medical Device Yes

Is Serialized Item No
Is Dangerous Good No
Shelf Life (Months) 60

## **PRODUCT DESCRIPTION**

AeroStep™ Protective Bunion Sleeve is an elastic protective sleeve. The inner gel pad provides relief and protection in the bunion area. It's designed to fit most types of footwear.

**DIRECTIONS:** Place over big toe, and position pad over the bunion. Wash with warm water and let dry. Use for more than 16 hours per day is not recommended.

TOPICAL USE ONLY - DO NOT PUT IN CONTACT WITH OPEN WOUNDS. DO NOT CLEAN WITH ALCOHOL. KEEP OUT OF REACH OF CHILDREN.

THIS PRODUCT IS NOT A SUBSTITUTE FOR PROFESSIONAL ATTENTION FROM A HEALTHCARE PRACTITIONER.

Location

DIABETICS AND PEOPLE WITH BLOOD CIRCULATION PROBLEMS SHOULD CONSULT WITH A PHYSICIAN BEFORE USE.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



## AEROSTEP PROTECTIVE BUNION SLEEVE 8.5 X 9CM BOX/1

ASTBS01