



ADDITIONAL INFORMATION

Weight	0.022 kg
Dimensions	14.7 × 7.6 × 2 cm
Carton Qty	585
Pricing Unit	Box/1
Min Order Qty	1
ARTG ID	426930
GMDN Code	36295
Medical Device Class	Class 1
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	9341394034228



PRODUCT DESCRIPTION

AeroStep™ Protective Bunion Sleeve is an elastic protective sleeve. The inner gel pad provides relief and protection in the bunion area. It's designed to fit most types of footwear.

DIRECTIONS: Place over big toe, and position pad over the bunion. Wash with warm water and let dry. Use for more than 16 hours per day is not recommended.

TOPICAL USE ONLY - DO NOT PUT IN CONTACT WITH OPEN WOUNDS. DO NOT CLEAN WITH ALCOHOL. KEEP OUT OF REACH OF CHILDREN.

THIS PRODUCT IS NOT A SUBSTITUTE FOR PROFESSIONAL ATTENTION FROM A HEALTHCARE PRACTITIONER.

DIABETICS AND PEOPLE WITH BLOOD CIRCULATION PROBLEMS SHOULD CONSULT WITH A PHYSICIAN BEFORE USE.



AEROSTEP PROTECTIVE BUNION SLEEVE 8.5 X 9CM BOX/1

ASTBS01