AERO Healthcare

DISPOSABLE GREEN PLASTIC KIDNEY DISH 200ML

ADDITIONAL INFORMATION



Weight	0.03 kg
Dimensions	16 × 10 × 4 cm
Carton Qty	100
Pricing Unit	Unit/1
Size	200ml
Min Order Qty	1
UNSPSC	42142529
Tax Schedule	Taxable
Is Medical Device	No
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

The AEROSUPPLIES stationary range covers a number of common stationary requirements applicable to the First Aid Industry.

Disposable Plastic Kidney Dish

Kidney dishes are a widely used product in surgical procedures for holding small instruments and more importantly for rinsing with sterile water. The dishes are also used to contain samples for laboratory analysis. The AeroSupplies Plastic Kidney dish is disposable.

We are able to source autoclavable kidney dishes if you have a specific request for reusable kidney dishes.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Plastic Kidney Dish

DISPOSABLE GREEN PLASTIC KIDNEY DISH 200ML



KD150D