



ADDITIONAL INFORMATION

| | |
|---------------------------|----------------|
| Weight | 0.03 kg |
| Dimensions | 16 × 10 × 4 cm |
| Carton Qty | 100 |
| Pricing Unit | Unit/1 |
| Size | 200ml |
| Min Order Qty | 1 |
| UNSPSC | 42142529 |
| Tax Schedule | Taxable |
| Is Serialized Item | No |
| Is Dangerous Good | No |

PRODUCT DESCRIPTION

The AEROSUPPLIES stationary range covers a number of common stationary requirements applicable to the First Aid Industry.

Disposable Plastic Kidney Dish

Kidney dishes are a widely used product in surgical procedures for holding small instruments and more importantly for rinsing with sterile water. The dishes are also used to contain samples for laboratory analysis. The AeroSupplies Plastic Kidney dish is disposable.

We are able to source autoclavable kidney dishes if you have a specific request for reusable kidney dishes.

Plastic Kidney Dish