



ADDITIONAL INFORMATION

Weight	0.08 kg
Dimensions	7.5 × 7.5 × 8.5 cm
Carton Qty	80
Pricing Unit	Pack/100
Size	7.5cm x 7.5cm
Min Order Qty	1
ARTG ID	236938
GMDN Code	11315
Medical Device Class	Class 1
UNSPSC	42311512
Tax Schedule	Taxable
Item Length	7.5
Item Width	7.5
Item Height	7.5
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	60

PRODUCT DESCRIPTION

The AEROSWAB™ Gauze Swab is highly absorbent cotton gauze swab that when soaked in fluid becomes an ideal item to clean and prepare wounds for dressing.

Cotton Gauze

- AEROSWAB™ cotton gauze has a balanced combination of softness, abrasive texture and absorbency for effectively



AEROSWAB NON-STERILE WHITE GAUZE SWAB 7.5 X 7.5CM PACK/100 (8PLY)

AGS75

cleaning wounds.

Multiple Sizes Available

- AEROSWAB™ gauze swabs are available in a number of different sizes to suit any application

8 ply Gauze Swab