

## AEROSWAB STERILE WHITE GAUZE SWAB 10 X 10CM (PACKS OF 3) BOX/25

**AGS103S** 







## **ADDITIONAL INFORMATION**

Weight 0.2 kg

**Dimensions**  $13 \times 9.5 \times 15.3 \text{ cm}$ 

Carton Qty 24

Pricing Unit Box/25

Size 10cm x 10cm

Min Order Qty 1

**ARTG ID** 159431 **GMDN Code** 13700

Medical Device Class 1s

**UNSPSC** 42311512

Tax Schedule Taxable

**Item Length** 10

**Item Width** 10

Item Measurement Unit cm

**Is Medical Device** Yes

**Is Serialized Item** No

No

Is Dangerous Good

## PRODUCT DESCRIPTION

AEROSWAB Sterile White Gauze Swab is a soft, 100% cotton gauze square designed for dependable absorbency when cleaning, blotting and cushioning minor wounds in first aid and clinical settings. The balanced texture supports effective wiping while being gentle on skin.

Sterile processing provides assurance for contact with clean wounds and procedural use. The open-weave cotton readily takes up saline and antiseptic solutions, making it suitable for wound preparation, light packing and protection of surrounding skin under secondary dressings and tapes.

Each swab measures 5 x 5 cm for easy handling and targeted coverage on hands, arms and smaller treatment sites. The crisp white colour helps visual checks during care, and the pliable cotton conforms smoothly to the body for comfortable application.

Trusted across healthcare, workplace kits and home first aid, AEROSWAB offers reliable performance and straightforward selection across common sizes, enabling consistent outcomes from initial clean-up through to dressing change support.



## AEROSWAB STERILE WHITE GAUZE SWAB 10 X 10CM (PACKS OF 3) BOX/25

**AGS103S** 

- 100% cotton gauze: soft, skin-friendly absorbency for cleaning and cushioning.
- Sterile: suitable for contact with clean wounds and procedural use.
- **High fluid uptake:** readily absorbs saline and exudate to support care.
- **5 x 5 cm format:** easy handling and focused coverage for small areas.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Gauze Swab - 3 Pcs