



ADDITIONAL INFORMATION

Weight	0.0064 kg
Dimensions	10.8 × 12.5 × 13 cm
Carton Qty	24
Pricing Unit	Box/25
Size	7.5cm x 7.5cm
Min Order Qty	1
ARTG ID	159431
GMDN Code	13700
Medical Device Class	Class 1s
UNSPSC	42311512
Tax Schedule	Taxable
Item Length	7.5
Item Width	7.5
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	60

PRODUCT DESCRIPTION

The AEROSWAB™ Gauze Swab is highly absorbent cotton gauze swab that when soaked in fluid becomes an ideal item to clean and prepare wounds for dressing.

Cotton Gauze

- AEROSWAB™ cotton gauze has a balanced combination of softness, abrasive texture and absorbency for effectively cleaning wounds.



AEROSWAB STERILE WHITE GAUZE SWAB 7.5 X 7.5CM (PACKS OF 5) BOX/25

AGS755S

Multiple Sizes Available

- AEROSWAB™ gauze swabs are available in a number of different sizes to suit any application.

Sterile

- AEROSWAB™ gauze swabs are also available in a sterile form for minimising infection risk when dressing clean wounds

Gauze Swab - 5 Pcs