

## AEROTAPE ELASTIC FABRIC ADHESIVE TAPE 5CM X

ATE50



## **ADDITIONAL INFORMATION**

Weight 0.05 kg

**Dimensions**  $5.8 \times 5.8 \times 5.8 \text{ cm}$ 

Carton Qty 36

Pricing Unit Each

Size 5cm x 2M

Min Order Qty 3

**ARTG ID** 236937 **GMDN Code** 34830

Medical Device Class 1

**UNSPSC** 42311552

Tax Schedule Taxable

Item Length 200

**Item Width** 5

Item Measurement Unit cm

**Is Medical Device** Yes

**Is Serialized Item** No

**Is Dangerous Good** No

Shelf Life (Months) 36

## **PRODUCT DESCRIPTION**

AeroTape Elastic Fabric Adhesive Fixation Tape is a versatile and comfortable tape designed for use in medical and sports settings. It is made from a flexible and breathable elastic fabric that conforms to the body, allowing for a full range of motion while providing support and stability. The adhesive is strong and durable, yet gentle on the skin, ensuring that the tape stays securely in place without causing irritation or pain. While being gentle on even the most sensitive skin, it can be easily removed without causing irritation or pain. This tape can be used to secure bandages, IVs, tubing, and other medical devices.

The 5cm x 5m tape provides a larger-sized Elastic Fabric Adhesive Tape for securely fixating items.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



## AEROTAPE ELASTIC FABRIC ADHESIVE TAPE 5CM X Healthcare

ATE50