



ADDITIONAL INFORMATION

Weight	0.24 kg
Dimensions	5.8 × 5.8 × 3.3 cm
Carton Qty	72
Pricing Unit	Each
Min Order Qty	6
ARTG ID	236937
GMDN Code	34830
Medical Device Class	Class 1
UNSPSC	42311552
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	9341394020870

PRODUCT DESCRIPTION

AeroTape Silk Fabric Adhesive Tape is a premium, hypoallergenic tape, ideal for use on sensitive skin. Made from a soft and breathable silk fabric that conforms to the body, providing comfort and flexibility. The adhesive is medical-grade, gentle, and skin-friendly and specially formulated to minimize irritation and pain during removal. This tape can be used to secure bandages, IVs, tubing, and other medical devices, as well as for wound care, and is perfect for hospitals, clinics, and home healthcare environments. The silk fabric's natural properties also make it ideal for use in situations where a discreet and aesthetically pleasing appearance is desired.

The 2.5cm x 5m tape provides a medium-sized Elastic Fabric Adhesive Tape for discreetly fixing items.