



## ADDITIONAL INFORMATION

<b>Weight</b>	0.135 kg
<b>Dimensions</b>	5.8 × 5.8 × 2 cm
<b>Carton Qty</b>	108
<b>Pricing Unit</b>	Each
<b>Size</b>	1.25cm x 5M
<b>Min Order Qty</b>	9
<b>ARTG ID</b>	236937
<b>GMDN Code</b>	34830
<b>Medical Device Class</b>	Class 1
<b>UNSPSC Code</b>	42311552
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	500
<b>Item Width</b>	1.25
<b>Item Measurement Unit</b>	cm
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Shelf Life (Months)</b>	36



## PRODUCT DESCRIPTION

AeroTape Waterproof Adhesive Tape is designed for reliable fixation in wet or humid conditions, helping keep dressings and medical devices securely in place when exposure to water or moisture is likely. This high tensile strength tape performs effectively in everyday first aid, clinical environments, and active settings where standard tapes may loosen or peel away. The waterproof construction supports firm fixation while maintaining flexibility so it can conform comfortably around limbs, joints, and awkward contours.



# AEROTAPE WATERPROOF ADHESIVE TAPE 1.25CM X 5M

ATW12

The tape uses a strong, medical-grade adhesive to provide dependable hold on the skin, even when the surface is exposed to sweat or splashes. This adhesive is hypoallergenic and latex free, making it suitable for users with latex sensitivity and helping to reduce the risk of skin reactions during extended wear. Despite its robust bond, the tape is designed to remove cleanly and with minimal discomfort, supporting patient comfort during dressing changes at home or in professional care.

Easy-tear properties allow the tape to be quickly torn by hand without the need for scissors, which saves time in fast-paced clinical workflows and everyday first aid. The narrow profile is well suited to securing gauze pads, lightweight dressings and small devices while minimising bulk on the skin. It is ideal for use in hospitals, clinics, sports and workplace first aid rooms, and for home healthcare where showers, handwashing and outdoor activities are part of daily life.

Across the range, AeroTape Waterproof Adhesive Tape supports secure fixation for bandages, IV lines, tubing and other devices that must stay in place despite contact with water. It offers a practical balance of durability, comfort and convenience, making it a useful addition to home first aid kits, clinical trolleys and emergency response bags.

- **Waterproof design:** maintains secure fixation in wet or humid environments.
- **Strong medical-grade adhesive:** provides reliable hold while remaining gentle on skin.
- **Latex free and hypoallergenic:** suitable for users with latex sensitivity.
- **Hand tearable construction:** enables quick, tool-free application and adjustment.
- **Flexible and conformable:** wraps comfortably around joints and contoured areas.

The **Size 1.25 cm x 5 m** option provides a narrow, easy-to-handle roll that is ideal for securing small dressings, gauze pads and tubing while keeping bulk to a minimum and allowing clothing to sit comfortably over the area.

This roll supplies **1.25 cm x 5 m** of waterproof adhesive tape, providing ample length for multiple small dressings while remaining compact enough to fit easily into first aid kits, medical bags and household supplies.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.