



ADDITIONAL INFORMATION

Weight	0.135 kg
Dimensions	5.8 × 5.8 × 2 cm
Carton Qty	108
Pricing Unit	Each
Size	1.25cm x 5M
Min Order Qty	9
ARTG ID	236937
GMDN Code	34830
Medical Device Class	Class 1
UNSPSC	42311552
Tax Schedule	Taxable
Item Length	500
Item Width	1.25
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	9341394020894

PRODUCT DESCRIPTION

AeroTape Waterproof Adhesive Fixation tape is a strong, durable tape that is designed specifically for use in wet or damp environments. Made from a waterproof, medical-grade adhesive designed to securely fixate bandages, IVs, tubing, and other medical devices to the skin. The tape is also easy to tear and apply, perfect for a wide range of applications. While being gentle on even the most sensitive skin, it can be easily removed without causing irritation or pain. AeroTape Waterproof is ideal for use in a variety of settings, including hospitals, clinics, and home healthcare environments.

The 1.25cm x 5m tape provides a narrower-sized Waterproof Adhesive Tape for securely fixing items.



AEROTAPE WATERPROOF ADHESIVE TAPE 1.25CM X 5M

ATW12