

AEROTAPE WATERPROOF ADHESIVE TAPE 2.5CM

ATW25





ADDITIONAL INFORMATION

Weight 0.04 kg

Dimensions $5.8 \times 5.8 \times 3.3 \text{ cm}$

Carton Qty 72
Pricing Unit Each

Min Order Qty 6

ARTG ID 236937
GMDN Code 34830
Medical Device Class 1

UNSPSC 42311552

Tax Schedule Taxable

Is Medical Device Yes
Is Serialized Item No
Is Dangerous Good No
Shelf Life (Months) 36

PRODUCT DESCRIPTION

AeroTape Waterproof Adhesive Fixation tape is a strong, durable tape that is designed specifically for use in wet or damp environments. Made from a waterproof, medical-grade adhesive designed to securely fixate bandages, IVs, tubing, and other medical devices to the skin. The tape is also easy to tear and apply, perfect for a wide range of applications. While being gentle on even the most sensitive skin, it can be easily removed without causing irritation or pain. AeroTape Waterproof is ideal for use in a variety of settings, including hospitals, clinics, and home healthcare environments.

The 2.5cm x 5m tape provides a medium-sized Waterproof Adhesive Tape for securely fixating items.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



AEROTAPE WATERPROOF ADHESIVE TAPE 2.5CM X 5M

ATW25