

AEROWASTE BIO-WASTE VOMIT BAG 1500ML

AW1514



ADDITIONAL INFORMATION

Weight 0.1 kg

Dimensions $5 \times 17.8 \times 21.5 \text{ cm}$

Carton Qty 48

Pricing Unit Bag/4

Size 1500ml

Min Order Qty 6

ARTG ID 236792
GMDN Code 36258
Medical Device Class 1

UNSPSC 47121707

Tax Schedule Taxable

Is Medical Device Yes
Is Serialized Item No
Is Dangerous Good No

PRODUCT DESCRIPTION

AEROWASTE products offer a convenient solution for collection of hazardous waste with the vomit bag for travel, clinical and hospital use.

Tough and Durable

The AEROWASTE Vomit Bag is manufactered from premium raw materials to ensure no spillage caused by structural failure of the bag

Twist Lock

The twist lock feature seales the bag, eliminating odour and risk of spillage

Measurement Gradation

The bag is clearly marked with measurement gradations to accurate measure the amount of bodily fluids lost

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



AEROWASTE BIO-WASTE VOMIT BAG 1500ML BAG/4

AW1514

Vomit bag with twist lock