

ADDITIONAL INFORMATION



Weight	7 kg
Dimensions	38 × 29 × 31 cm
Carton Qty	1
Pricing Unit	Box/2000
Size	200mm x 100mm
Min Order Qty	1
ARTG ID	277597
GMDN Code	61694
Medical Device Class	Class 1
UNSPSC	42141503
Tax Schedule	Taxable
Item Length	20
Item Width	10
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

The AEROWIPE cleansing wipes are a moist tissue containing an aqueous solution with 1% cetrimide BP W/V and 0.02% Chlorhexidine BP W/V.

These are used to soothe and clean cuts and grazes, helping to maximise the prevention of infection.

This skin cleansing product is alcohol free, making it patient friendly with the absence of fumes and the stinging sensation caused by conventional alcohol wipes.



AEROWIPE ALCOHOL-FREE CLEANSING WIPES BOX/2000

AW8000

Alcohol-Free Cleansing Wipes