

AEROWIPE ALCOHOL-FREE CLEANSING WIPES ENV/10

AW8110





ADDITIONAL INFORMATION

Weight 0.03889 kg

Dimensions $7.5 \times 8 \times 0.5 \text{ cm}$

Carton Qty 200

Pricing Unit Env/10

Size 200mm x 100mm

Min Order Qty 200

ARTG ID 277597

GMDN Code 61694

Medical Device Class 1

UNSPSC 42141503

Tax Schedule Taxable

Item Length 20

Item Width 10

Item Measurement Unit cm

Is Medical Device Yes

Is Serialized Item No

Is Dangerous Good No

PRODUCT DESCRIPTION

Introducing the AEROWIPE Alcohol-Free Cleansing Wipes, an essential addition to any first aid kit. These wipes are expertly designed to provide a gentle yet effective solution for maintaining skin hygiene and cleanliness. Made from high-quality absorbent fabric, each wipe is pre-soaked in a carefully formulated aqueous solution containing 1% cetrimide BP W/V and 0.02% Chlorhexidine BP W/V. This unique blend ensures thorough cleansing while being kind to the skin, making them suitable for a wide range of applications.

One of the key advantages of these wipes is their alcohol-free composition, which eliminates the discomfort of stinging sensations and harsh fumes often associated with alcohol-based products. This makes them particularly suitable for individuals with sensitive skin.

Designed as non-sterile, single-use devices, these wipes are intended for cleaning localized body surfaces before medical procedures or device applications. They are also perfect for everyday skin hygiene, providing a convenient and effective way to keep skin clean and healthy. Available over-the-counter, AEROWIPE Alcohol-Free Cleansing Wipes are a practical choice for anyone seeking a dependable and gentle cleansing solution.



AEROWIPE ALCOHOL-FREE CLEANSING WIPES ENV/10

AW8110

This	product is not a s	substitute for	direct medical	attention fr	om a healthca	are practitioner	. Always read t	the label ar	nd follow the
dire	ctions for use.					·	-		

Alcohol-Free Cleansing Wipes