



ADDITIONAL INFORMATION

Weight	0.03889 kg
Dimensions	7.5 × 8 × 0.5 cm
Carton Qty	200
Pricing Unit	Env/10
Size	200mm x 100mm
Min Order Qty	200
ARTG ID	277597
GMDN Code	61694
Medical Device Class	Class 1
UNSPSC Code	42141503
Tax Schedule	Taxable
Item Length	20
Item Width	10
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

AEROWIPE Alcohol-Free Cleansing Wipes are designed for gentle, effective skin hygiene in first aid kits, workplaces, clinical environments and everyday routines. Each single-use wipe is made from strong non-woven fabric and pre-soaked with an aqueous solution containing **1% cetrimide BP W/V** and **0.02% chlorhexidine BP W/V**. This formulation supports dependable cleansing of localised areas of skin before applying dressings or devices, and as part of routine first aid care.

The alcohol-free formulation helps avoid the sting, strong odour and fumes often associated with alcohol-based wipes, making it a

comfortable option for frequent use and for people who prefer a gentler wipe. The non-woven carrier is durable and low-lint, helping the wipe glide smoothly and maintain a tidy clean-up. The wipe size is generous for efficient coverage of small areas without needing multiple wipes for a single task.

Individually wrapped sachets help keep each wipe hygienic until needed, making it easy to carry in a bag, glove box or kit, and simple to restock first aid supplies. AEROWIPE Alcohol-Free Cleansing Wipes are non-sterile, intended for external use, and suited to quick skin cleansing when soap and water are not readily available.

- **Alcohol-free comfort:** helps avoid sting and strong odour.
- **Dual formulation:** 1% cetrimide and 0.02% chlorhexidine for effective cleansing.
- **Strong non-woven wipe:** durable, low-lint carrier for controlled wipe-downs.
- **Individually wrapped:** hygienic, portable and convenient for kits and travel.
- **Single-use and non-sterile:** practical format for routine first aid hygiene.

Quantity: Bag/10 individually wrapped wipes is ideal for keeping a number of AEROWIPE Alcohol-Free Cleansing Wipes handy in a travel bag, glovebox, or compact first aid kit, while still providing the same alcohol-free skin cleansing performance.

Supplied in a envelope of **10** individually wrapped wipes measuring **200 mm x 100 mm**, providing broad coverage for efficient cleansing and straightforward first aid kit replenishment.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Alcohol-Free Cleansing Wipes