

## AEROWIPE 70% ISOPROPYL ALCOHOL SWAB 3 X 3CM BOX/100

AW7100





## **ADDITIONAL INFORMATION**

Weight 0.1 kg

**Dimensions**  $5.4 \times 10.2 \times 5.4 \text{ cm}$ 

Carton Qty 100

Pricing Unit Box/100

**Size** 30 x 30mm

Min Order Qty 1

ARTG ID 277597
GMDN Code 61694
Medical Device Class 1

**UNSPSC** 42141503

Tax Schedule Non-Taxable

Item Length30Item Width3

**Item Measurement Unit** cm

Is Medical Device Yes
Is Serialized Item No
Is Dangerous Good No

60

Shelf Life (Months)

## **PRODUCT DESCRIPTION**

AeroWipe Skin prep wipes contain 70% Isopropyl Alcohol which is an effective antibacterial formula for cleansing the skin prior to injection

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Alcohol Swabs (2ply)



## AEROWIPE 70% ISOPROPYL ALCOHOL SWAB 3 X 3CM BOX/100

AW7100