



ADDITIONAL INFORMATION

Weight	0.17 kg
Dimensions	7 × 10.2 × 6 cm
Carton Qty	100
Pricing Unit	Box/100
Size	60mm x 33mm
Min Order Qty	1
ARTG ID	277597
GMDN Code	61694
Medical Device Class	Class 1
UNSPSC Code	42141503
Tax Schedule	Taxable
Item Length	6
Item Width	3.3
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

AEROWIPE 10% Povidone Iodine Swabs deliver dependable pre-procedure skin cleansing for homes, clinics and workplaces. The absorbent fabric pad is pre-saturated with povidone iodine to create a uniform brown film that clearly shows coverage, helping users prepare a localised site before injections, dressings, vaccinations and minor clinical tasks. The alcohol-free composition makes it a comfortable choice for intact skin while still providing trusted cleansing action.

The wipe's 2-ply construction combines strength and absorbency so solution spreads evenly across both flat and contoured areas



AEROWIPE 10% POVIDONE IODINE SWABS 60 X 33MM BOX/100

AWP16100

without excessive pooling. A compact 3 cm x 6 cm format sits securely between the fingers for precise control around shoulders, arms and legs, supporting quick, consistent preparation in busy care settings. Each swab is non-sterile and designed for single use to support good hygiene practice between patients and steps in a procedure.

Individually wrapped sachets keep each pad ready to use in first aid kits, treatment rooms, vaccination stations and workplace medical cabinets. The visible iodine tint helps users verify coverage at a glance and promotes a repeatable technique across teams and locations. From day-to-day first aid to routine clinic workflows, this swab supports efficient, confident site preparation that fits neatly into standard operating procedures.

- **10% povidone iodine:** reliable skin cleansing.
- **Alcohol free:** comfortable on intact skin.
- **2-ply pad:** strong, absorbent and even coverage.
- **3 cm x 6 cm size:** compact control for precise application.
- **Single use, non-sterile:** supports good hygiene practice.

Supplied as a box of **100** individually wrapped, alcohol-free swabs saturated with **10% povidone iodine** to standardise pre-procedure cleansing; the compact wipe format supports quick, even coverage.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Povidone Iodine Swabs (2 ply)