



ADDITIONAL INFORMATION

Weight	0.08 kg
Dimensions	14.7 × 3.5 × 16.5 cm
Carton Qty	20
Pricing Unit	Box/10
Min Order Qty	1
ARTG ID	406570
GMDN Code	43186
Medical Device Class	Class 2b
Tax Schedule	Non-Taxable
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

AeroWound Alginate Dressings provide an advanced wound care solution, designed to provide a moist healing environment. Naturally derived from seaweed, AeroWound is an ideal choice for various wound types. On contact with wound exudate, the dressing transforms into a soft gel, maintaining a moist healing environment which is conducive to wound healing, and assists in pain reduction.

<https://aerohealthcareonline.com/wp-content/uploads/2023/04/AlginateDressing-Final.mp4>

Features:

Optimized Healing Environment: Our dressings ensure a moist wound environment, speeding up the natural healing process.

High Absorbency: Specially crafted to absorb wound exudate efficiently, minimising the risk of maceration.

Natural Seaweed Derivative: Made from natural seaweed, providing bio-degradable, eco-friendly care.

Soft & Conforming: Flexible design adheres comfortably to the wound contours, ensuring patient comfort.



AEROWOUND ALGINATE DRESSING 10 X 10CM BOX/10

AWA010

Easy Removal: Designed to gel up upon contact with exudate for painless removal, reducing patient trauma.

Safety Considerations: For external use only. Always read the label and follow the directions for use.

Alginate Dressing 5 x 5cm