



## ADDITIONAL INFORMATION

<b>Weight</b>	0.08 kg
<b>Dimensions</b>	14.7 × 3.5 × 16.5 cm
<b>Carton Qty</b>	20
<b>Pricing Unit</b>	Box/10
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	406570
<b>GMDN Code</b>	43186
<b>Medical Device Class</b>	Class 2b
<b>Tax Schedule</b>	Non-Taxable
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No

## PRODUCT DESCRIPTION

AeroWound Alginate Dressings provide an advanced wound care solution, designed to provide a moist healing environment. Naturally derived from seaweed, AeroWound is an ideal choice for various wound types. On contact with wound exudate, the dressing transforms into a soft gel, maintaining a moist healing environment which is conducive to wound healing, and assists in pain reduction.

<https://aerohealthcareonline.com/wp-content/uploads/2023/04/AlginateDressing-Final.mp4>

### Features:

**Optimized Healing Environment:** Our dressings ensure a moist wound environment, speeding up the natural healing process.

**High Absorbency:** Specially crafted to absorb wound exudate efficiently, minimising the risk of maceration.

**Natural Seaweed Derivative:** Made from natural seaweed, providing bio-degradable, eco-friendly care.

**Soft & Conforming:** Flexible design adheres comfortably to the wound contours, ensuring patient comfort.



# AEROWOUND ALGINATE DRESSING 10 X 10CM BOX/10

AWA010

**Easy Removal:** Designed to gel up upon contact with exudate for painless removal, reducing patient trauma.

**Safety Considerations:** For external use only. Always read the label and follow the directions for use.

Alginate Dressing 5 x 5cm